



Muse Developer Resources

Search this site

Navigation

Home
SDK Overview
Download & Install
Getting Started
LibMuse
MuseIO
MuseLab
MusePlayer
MuseIO Receiver
Muse Hardware
Muse Data Files
Muse Communication Protocol
Multi-Muse Setup
Developer FAQ
Intro to BCI and EEG
Release Notes
Forums
Mailing List
Support

MusePlayer

Subpage Listing

[MATLAB Output File](#)

Overview

[MusePlayer](#) is a utility for recording, replaying, rerouting, and converting EEG and accelerometer data from Interaxon Muse EEG devices. It can save to and convert between the native Muse datatype (.muse), Matlab (HDF5), CSV, and OSC replay formats.

It is an open source project. All the source code and examples are available on [Bitbucket](#).

Currently muse-player supports the following inputs and outputs:

Supported inputs:

- OSC network stream
- OSC-replay file format
- Muse file format v1
- Muse file format v2

Supported outputs:

- MATLAB ([click here for a description of the output](#))
- CSV
- OSC network stream
- OSC-replay file format
- Muse file format v2
- Print to screen

How to run it

Simply type "muse-player" followed by your desired input and output arguments. See below for a description of these arguments. You can also view this information by running muse-player with no arguments - it will be printed to your screen.

optional arguments:

`-h, --help` show this help message and exit

`-v, --verbose` Print status messages to stdout

`-q, --as-fast-as-possible` Replay input as fast as possible instead of using original timing.

`-j, --jump-data-gaps` Replay input by omitting any data gaps larger than 1 second.

`-n, --no--time--data` Replay input by omitting output of current timing info.

`-i FILTER_DATA [FILTER_DATA ...], --filter FILTER_DATA [FILTER_DATA ...]` Filter data by path. e.g.
`-i /muse/dsp/elements/alpha`
`/muse/eeg`

Input options:

Only one type of input can be specified, but can

be multiple files of the same type:

```
-l [INPUT_OSC_PORT], --input-osc-port  
[INPUT_OSC_PORT]
```

Listen for OSC messages on
this port (default:

tcp:5000).

```
-f INPUT_MUSE_FILES [INPUT_MUSE_FILES ...],  
--input-muse-files INPUT_MUSE_FILES  
[INPUT_MUSE_FILES ...]
```

Input from Muse file
format.

```
-o INPUT_OSCREPLAY_FILES [INPUT_OSCREPLAY_FILES  
...], --input-oscreplay-files INPUT_OSCREPLAY_FILES  
[INPUT_OSCREPLAY_FILES ...]
```

Input from OSC-replay
files.

Output options:

One or more outputs can be specified:

```
-s [OUTPUT_OSC_URL], --output-osc-url  
[OUTPUT_OSC_URL]
```

Output OSC messages to
HOST:PORT (default:

osc.tcp://localhost:5000)

```
-F FILE, --output-muse-file FILE
```

Output to a Muse file

```
-M FILE, --output-matlab-file FILE
```

Output to a Matlab file

`-O FILE, --output-oscreplay-file FILE`

Output to an OSC-replay

file

`-C FILE, --output-csv-file FILE`

Output to an CSV file

`-D, --output-screen-dump`

Output to the screen

directly

Examples:

```
muse-player -f my_eeg_recording.muse -s  
osc.tcp://localhost:7887
```

This will read in the file
"my_eeg_recording.muse" and send those messages as
OSC to port 7887.

```
muse-player -l 5555 -M matlab.mat -s 5000
```

This will receive OSC messages on port
5555, save them to file, and rebroadcast them to
port 5000.



Subpages (1): [MATLAB Output File](#)

You do not have permission to add comments.

Comments