

 Search this site

## Navigation

Home  
SDK Overview  
Download & Install  
Getting Started  
LibMuse  
MuseIO  
MuseLab  
**MusePlayer**  
MuseIO Receiver  
Muse Hardware  
Muse Data Files  
Muse Communication Protocol  
Multi-Muse Setup  
Developer FAQ  
Intro to BCI and EEG  
Release Notes  
Forums  
Mailing List  
Support

## MusePlayer

### Subpage Listing

[MATLAB Output File](#)

## Overview

[MusePlayer](#) is a utility for recording, replaying, rerouting, and converting EEG and accelerometer data from Interaxon Muse EEG devices. It can save to and convert between the native Muse datatype (.muse), Matlab (HDF5), CSV, and OSC replay formats.

It is an open source project. All the source code and examples are available on [Bitbucket](#).

Currently muse-player supports the following inputs and outputs:

Supported inputs:

- OSC network stream
- OSC-replay file format
- Muse file format v1
- Muse file format v2

Supported outputs:

- MATLAB ([click here for a description of the output](#))
- CSV
- OSC network stream
- OSC-replay file format
- Muse file format v2
- Print to screen

## How to run it

Simply type "muse-player" followed by your desired input and output

arguments. See below for a description of these arguments. You can also view this information by running `muse-player` with no arguments - it will be printed to your screen.

optional arguments:

<code>-h, --help</code>	show this help message and exit
<code>-v, --verbose</code>	Print status messages to stdout
<code>-q, --as-fast-as-possible</code>	Replay input as fast as possible instead of using original timing.
<code>-j, --jump-data-gaps</code>	Replay input by omitting any data gaps larger than 1 second.
<code>-n, --no--time--data</code>	Replay input by omitting output of current timing info.
<code>-i FILTER_DATA [FILTER_DATA ...], --filter FILTER_DATA [FILTER_DATA ...]</code>	Filter data by path. e.g. <code>-i /muse/dsp/elements/alpha</code> <code>/muse/eeg</code>

Input options:

Only one type of input can be specified, but can be multiple files of the same type:

<code>-l [INPUT_OSC_PORT], --input-osc-port [INPUT_OSC_PORT]</code>	Listen for OSC messages on this port (default: <code>tcp:5000</code> ).
---	---

```
-f INPUT_MUSE_FILES [INPUT_MUSE_FILES ...], --input-  
muse-files INPUT_MUSE_FILES [INPUT_MUSE_FILES ...]
```

Input from Muse file format.

```
-o INPUT_OSCREPLAY_FILES [INPUT_OSCREPLAY_FILES ...],  
--input-oscreplay-files INPUT_OSCREPLAY_FILES  
[INPUT_OSCREPLAY_FILES ...]
```

Input from OSC-replay files.

#### Output options:

One or more outputs can be specified:

```
-s [OUTPUT_OSC_URL], --output-osc-url [OUTPUT_OSC_URL]
```

Output OSC messages to HOST:PORT

(default:

```
osc.tcp://localhost:5000)
```

```
-F FILE, --output-muse-file FILE
```

Output to a Muse file

```
-M FILE, --output-matlab-file FILE
```

Output to a Matlab file

```
-O FILE, --output-oscreplay-file FILE
```

Output to an OSC-replay file

```
-C FILE, --output-csv-file FILE
```

Output to an CSV file

```
-D, --output-screen-dump
```

Output to the screen directly

#### Examples:

```
muse-player -f my_eeg_recording.muse -s  
osc.tcp://localhost:7887
```

```
This will read in the file  
"my_eeg_recording.muse" and send those messages as OSC to  
port 7887.
```

```
muse-player -l 5555 -M matlab.mat -s 5000
```

```
This will receive OSC messages on port 5555, save  
them to file, and rebroadcast them to port 5000.
```



---

Subpages (1): [MATLAB Output File](#)

## Comments

You do not have permission to add comments.